



insight iS3
Ultralight™ Serie

Light Reduces your operations burden

Omnipotent Satisfy your clinical demands



Insighters Medical

Add: Floor 13th, HTM Building, Ganli Road, Longgang District, Shenzhen, China
Tel: +86-755-89313133 ext 8128
Email: intl-market@insighters.cn
Website: www.insighters.cn

Features&Benefits

Light Design

The weight is only 250g, which reduces the burden on operators when holding the handle for a long time.

HD imaging

Full digital micro-imaging, no built-in optical fiber, high resolution imaging.

Wireless Connection

Support wireless connection with i-Workstation, release doctor's operation from cable limit .


Complete sizes portfolio for your needs

From iS3-C (5.8mm/3.0mm, 5.2mm/2.6mm, 4.5mm/2.0mm, 3.5mm/1.2mm) to iS3-F(2.4mm), the multiple diameter ranges fulfill your demands in various diagnostic and treatment applications.




Precise Navigation

Wide range up/down bending angle (iS3-C: $\geq 130^\circ/130^\circ$, iS3-F $\geq 150^\circ/150^\circ$, please contact insightfulers or our authorized distributors if higher bending angles needed), plus left/right direction adjustment function, ensure you to navigate inside patient's airway precisely and confidently.




Three levels of brightness adjustable

Meet the clinical needs of light source brightness in different scenarios, and assist in the diagnosis and treatment of diseases.



One-click freeze

Allow observe the clarity of the frozen picture in advance and obtain a guaranteed image, without affecting the real-time endoscopic inspection.



Rotation control

Just one button for quick and precise left/right direction adjustment.



Insertion Portion

Detailed front view of the distal end.

