

ORION™ II CT CVC

CVC Access with
Power Injectability



HEALTH LINE
MEDICAL PRODUCTS

Why choose a CT CVC?

Using an indwelling central line indicated for pressure injection of contrast media, negates the need for another type of venous access. This streamlines the scan process and saves crucial time for you and your patients.

The Orion™ II CT CVC catheters are power injectable and intended for short term infusion therapy. Available in a variety of sizes and lengths, clinicians have the options necessary to meet the patient's needs for central vein access.

Indications

- Blood Letting
- Fluid Monitoring
- Multiple Infusions
- Poor Peripheral Access
- Injection of Contrast Media
- Necessary Extended Therapy

Features and Benefits

Available in a Variety of Lumen Configurations - Offered in a single, double, and triple lumen.

Soft Flexible Tip - Reduces vessel irritation during insertion.

Reverse Taper Design - Strengthens the external portion of the catheter.

Simple Grip Luers - Easy on gloves and color coded for multiple treatment options.

One Piece Cannula Construction - No bonds or welds helps to ensure catheter will not break.

CVP Monitoring - ORION™ II Power CVCs are indicated for central vein pressure monitoring.

Polyurethan Catheter Construction - Provides radiopacity, durability, and strength.

Maximum Injection Rates - Allows for the injection of contrast media at rates up to 5 ml/sec for CECT scans.



Simple Grip Luers

are easy on gloves and maintain strength and crack resistance



Low Profile Hub
allows for suture or sutureless securement

Easy, Legible Identification
on clamps and extensions. Indicating power injectability and maximum pressure

Orion™ II CT CVC Kit Configurations

Catalogue No.	CT CVC Description	Guidewire	Box Qty.	Carton Qty.
A33-05113	16Ga x 13cm Single Lumen CT CVC	Stainless	5	5
A33-05215	5f x 15cm Dual Lumen CT CVC	Stainless	5	5
A33-05220	5f x 20cm Dual Lumen CT CVC	Stainless	5	5
A33-07215	7f x 15cm Dual Lumen CT CVC	Stainless	5	5
A33-07220	7f x 20cm Dual Lumen CT CVC	Stainless	5	5
A33-07315	7f x 15cm Triple Lumen CT CVC	Stainless	5	5
A33-07320	7f x 20cm Triple Lumen CT CVC	Stainless	5	5

Standard Tray Components

Component	Quantity	Component	Quantity
CT CVC Polyurethane Catheter	1	Syringe (10ml)	1
ICU MicroClave	1, 2, or 3	Raulerson Syringe (5ml)	1
.035" Guidewire	1	IFU	1
18Ga x 7cm Introducer Needle	1	Removable Suture Wing	2
8Fr x 10cm Dilator	1	Scalpel	1

For ordering information or for questions about Orion II, please contact your Health Line Sales Consultant at (877) 847-4542 or customerservice@hlic.net and visit us at www.hlic.net.



