

HEALTH LINE


CT Midline Catheters with
MST Insertion Kit



HEALTH LINE
MEDICAL PRODUCTS

The Health Line CT Midline Catheter is intended for short term (less than 30 days) peripheral access to the venous system for the purpose of intravenous therapy, medicines, and blood products. A single venipuncture can meet infusion therapy requirements for patients needing 5-7 days of therapy, thus eliminating multiple needle sticks and enhancing patient comfort. Additionally, the Health Line can reduce vein inflammation better than short peripheral catheters and does not require an X-ray after placement.

The Health Line is the ideal choice for clinicians and patients looking for reliable short-term access and is available in single and dual lumen configurations.



Features and Benefits

Depth Markings - Allow the catheter to be measured and trimmed accurately.

Kit Components - CT injectable Midline Catheter includes MST insertion components.

Polyurethane Catheter Construction - Provides patient with radiopacity, durability, and strength while still allowing for catheter to be soft inside the vessel.

Compatibility - Catheter is MRI compatible.

Placement - X-ray not required after placement.

Health Line CT Midline Catheters

Catalogue No.	PICC Description	Guide Wire	Introducer	Box Qty.	Carton Qty.
HEALTH LINE CT MIDLINE 20CM CATHETERS – NURSES KIT					
A120-03120-N	3f x 20cm Single Lumen CT Midline	40cm Nitinol GW w/ Stainless Steel Tip	3.5f x 05cm PTFE Tearaway	5	20
A120-04120-N	4f x 20cm Single Lumen CT Midline	40cm Nitinol GW w/ Stainless Steel Tip	4.5f x 05cm PTFE Tearaway	5	20
A120-04220-N	4f x 20cm Dual Lumen CT Midline	40cm Nitinol GW w/ Stainless Steel Tip	4.5f x 05cm PTFE Tearaway	5	20
A120-05120-N	5f x 20cm Single Lumen CT Midline	40cm Nitinol GW w/ Stainless Steel Tip	5.5f x 05cm PTFE Tearaway	5	20
A120-05220-N	5f x 20cm Dual Lumen CT Midline	40cm Nitinol GW w/ Stainless Steel Tip	5.5f x 05cm PTFE Tearaway	5	20



Standard Kit Components

Component	Quantity
Health Line CT Midline Catheter	1
Injection Cap	1, or 2
Twisted Wire Stylet	1
.030 I.D. Adaptor w/Sideport	1
21GA X 7CM Safety Introducer Needle	1
PTFE Tearaway Introducer	1
Nitinol Straight Mandrel GW w/ Stainless Steel Tip	1
10cc Disposable Luer Lock Syringe	1
Non-Safety Scalpel #11	1
24" Paper Tape Measure	1
Stylet TAG	1
Patient Info Guide	1
Chart Sticker	1
IFU	1
Securment Device	1
Midline Identification Band	1

For ordering information or for questions about the Health Line CT Midline, please contact your Health Line Sales Consultant at (877) 847-4542 and visit us at www.hlic.net.





Health Line International Corp.
260 North Ace Yeager Court, Unit D
Salt Lake City, Utah 84116

Phone: (877) 847-4542
Fax: (855) 228-1336
www.hlic.net

ISO13485
Certified Quality System
MKT-100-09 Rev. 01