

# ORION™

## Central Vein Catheters and Kits



Health Line International Corp.  
260 North Ace Yeager Court, Unit D  
Salt Lake City, Utah 84116

Phone: (877) 847-4542  
Fax: (855) 228-1336  
[www.hlic.net](http://www.hlic.net)

**ISO13485**  
**Certified Quality System**  
**MKT-200-02 Rev. 01**



# Orion™ Short Term Infusion Catheters and Kits

	Kit	Catalogue Number	Description	Guide Wire	Introducer Needle
SINGLE LUMEN CENTRAL VEIN CATHETER (CVC)	14Ga ORION™ Peripherally Inserted Central Catheters: Kit	A32-06110	14Ga x 10cm Single Lumen Orion™ CVC	Stainless	18Ga
		A32-06116	14Ga x 16cm Single Lumen Orion™ CVC	Stainless	18Ga
		A32-06120	14Ga x 20cm Single Lumen Orion™ CVC	Stainless	18Ga
		A32-06130	14Ga x 30cm Single Lumen Orion™ CVC	Stainless	18Ga
	16Ga ORION™ Peripherally Inserted Central Catheters: Kit	A32-05110	16Ga x 10cm Single Lumen Orion™ CVC	Stainless	18Ga
		A32-05116	16Ga x 16cm Single Lumen Orion™ CVC	Stainless	18Ga
		A32-05120	16Ga x 20cm Single Lumen Orion™ CVC	Stainless	18Ga
		A32-05130	16Ga x 30cm Single Lumen Orion™ CVC	Stainless	18Ga
	18Ga ORION™ Peripherally Inserted Central Catheters: Kit	A32-04110	18Ga x 10cm Single Lumen Orion™ CVC	Stainless	19Ga
		A32-04116	18Ga x 16cm Single Lumen Orion™ CVC	Stainless	19Ga
		A32-04120	18Ga x 20cm Single Lumen Orion™ CVC	Stainless	19Ga
		A32-04130	18Ga x 30cm Single Lumen Orion™ CVC	Stainless	19Ga
DOUBLE LUMEN CENTRAL VEIN CATHETER (CVC)	4fr ORION™ Peripherally Inserted Central Catheters: Kit	A32-04205	4f x 05cm Double Lumen Orion™ CVC	Stainless	20Ga
		A32-04208	4f x 08cm Double Lumen Orion™ CVC	Stainless	20Ga
		A32-04213	4f x 13cm Double Lumen Orion™ CVC	Stainless	20Ga
		A32-04230	4f x 30cm Double Lumen Orion™ CVC	Stainless	20Ga
	5fr ORION™ Peripherally Inserted Central Catheters: Kit	A32-05206	5f x 06cm Double Lumen Orion™ CVC	Stainless	19Ga
		A32-05208	5f x 08cm Double Lumen Orion™ CVC	Stainless	19Ga
		A32-05213	5f x 13cm Double Lumen Orion™ CVC	Stainless	19Ga
		7fr ORION™ Peripherally Inserted Central Catheters: Kit	A32-07215	7f x 15cm Double Lumen Orion™ CVC	Stainless
A32-07220	7f x 20cm Double Lumen Orion™ CVC		Stainless	18GA	
A32-07230	7f x 30cm Double Lumen Orion™ CVC		Stainless	18GA	
TRIPLE LUMEN CENTRAL VEIN CATHETER (CVC)	5.5fr ORION™ Peripherally Inserted Central Catheters: Kit		A32-05308	5.5f x 08cm Triple Lumen Orion™ CVC	Stainless
		A32-05313	5.5f x 13cm Triple Lumen Orion™ CVC	Stainless	19Ga
		A32-05330	5.5f x 30cm Triple Lumen Orion™ CVC	Stainless	19Ga
	7fr ORION™ Peripherally Inserted Central Catheters: Kit	A32-07315	7f x 15cm Triple Lumen Orion™ CVC	Stainless	18Ga
		A32-07320	7f x 20cm Triple Lumen Orion™ CVC	Stainless	18Ga
		A32-07230	7f x 30cm Triple Lumen Orion™ CVC	Stainless	18Ga


## Features and Benefits

**Thermal Reactive Polyurethane Material** – Softens after insertion for improved patient comfort.

**Simple Grip Luers** – Easy on gloves and color coded for multiple treatment options.

**Large Side Holes** – Improves flow rates and infusion times.

The ORION™ CVC is indicated for short term (less than 30 days) access to the central venous system for intravenous administration. They are available in a variety of sizes and lengths so that clinicians have the options necessary to meet the patient's needs.



ALL KITS CONTAIN THE FOLLOWING:

Component	Quantity	Component	Quantity
Orion™ Central Vein Catheter	1	Hypodermic Needle	1
Injection Cap	1, 2, or 3	10ml Syringe	1
Guidewire	1	5ml Raulerson Syringe	1
Introducer Needle	1	Movable Wing	1
Dilator	1	CVC IFU	1
#11 Scalpel	1		

For ordering information or for questions about the Orion™, please contact your Health Line Sales Consultant at (877) 847-4542 and visit us at [www.hlic.net](http://www.hlic.net).

