

# HEALTH LINE CT CVC

CVC Access with  
Power Injectability



## Why choose a CT CVC?

Using an indwelling central line indicated for pressure injection of contrast media, negates the need for another type of venous access. This streamlines the scan process and saves crucial time for you and your patients.

The Health Line CT CVC catheters are power injectable and intended for short term infusion therapy. Available in a variety of sizes and lengths, clinicians have the options necessary to meet the patient's needs for central vein access.

## Indications

Blood Letting  
Fluid Monitoring  
Multiple Infusions  
Poor Peripheral Access  
Injection of Contrast  
Media Necessary  
Extended Therapy

## Features and Benefits

**Available in a Variety of Lumen Configurations** – Offered in a single, double, and triple lumen.

**Soft Flexible Tip** - Reduces vessel irritation during insertion.

**Reverse Taper Design** - Strengthens the external portion of the catheter.

**Simple Grip Luers** – Easy on gloves and color coded for multiple treatment options.

**One Piece Cannula Construction** – No bonds or welds helps to ensure catheter will not break.

**CVP Monitoring** - ORION™ II Power CVCs are indicated for central vein pressure monitoring.

**Polyurethan Catheter Construction** - Provides radiopacity, durability, and strength.

**Maximum Injection Rates** - Allows for the injection of contrast media at rates up to 5 ml/sec for CECT scans.



**Simple Grip Luers**

are easy on gloves and maintain strength and crack resistance



**Low Profile Hub** allows for suture or sutureless securement

**Easy, Legible Identification**

on clamps and extensions. Indicating power injectability and maximum pressure

## Health Line CT CVC Kit Configurations

Catalogue No.	CT CVC Description	Guidewire	Box Qty	Carton Qty.
A33-05113	16 Ga x 13cm Single Lumen CT CVC	Stainless	5	20
A33-05215	5f x 15cm Dual Lumen CT CVC	Stainless	5	20
A33-05220	5f x 20cm Dual Lumen CT CVC	Stainless	5	20
A33-07215	7f x 15cm Dual Lumen CT CVC	Stainless	5	20
A33-07220	7f x 20cm Dual Lumen CT CVC	Stainless	5	20
A33-07315	7f x 15cm Triple Lumen CT CVC	Stainless	5	20
A33-07320	7f x 20cm Triple Lumen CT CVC	Stainless	5	20

## Standard Tray Components

Component	Quantity	Component	Quantity
CT CVC Polyurethane Catheter	1	Syringe (10ml)	1
ICU MicroClave	1, 2, or 3	Raulerson Syringe (5ml)	1
.035" Guidewire	1	IFU	1
18Ga x 7cm Introducer Needle	1	Removable Suture Wing	2
8Fr x 10cm Dilator	1	Scalpel	1

For ordering information or for questions about Health Line CT CVC, please contact your Health Line Sales Consultant at (877) 847-4542 or customerservice@hlic.net and visit us at www.hlic.net.





Health Line International Corp.  
260 North Ace Yeager Court, Unit D  
Salt Lake City, Utah 84116

Phone: (877) 847-4542  
Fax: (855) 228-1336  
[www.hlic.net](http://www.hlic.net)

**ISO13485**  
**Certified Quality System**  
**MKT-100-09 Rev. 01**