

THEATRE DRAPES & SURGICAL ACCESSORIES

DJ

ABOUT US

P3 Medical Limited is a British Company specialising in the design, manufacture and marketing of a range of high quality medical devices primarily for use within the operating theatre.

Our mission is to provide cost-effective solutions to enhance hospital efficiency, patient care and overall safety in operating theatres in the UK and abroad.

Within our two UK manufacturing sites in Bristol (Headquarters) and Preston, our facilities include 4 modern cleanrooms. We are proud to supply both healthcare organisations and international distributors with British products, manufactured to the highest standards. We currently export to more than 50 countries worldwide.

We are also an established UK distributor for a number of well-known international brands. These products complement our own product ranges and enable us to offer a diverse array of high quality products throughout the healthcare industry.

P3 Medical, established in 2001, is a privately owned Company with a wealth of experience in providing excellent service both directly to hospitals and via UK and international supply chains.

THEATRE DRAPES & SURGICAL ACCESSORIES

WELCOME D

OUR PRODUCTS

All our devices comply with the requirements of the Medical Device Directive and we operate a quality management system to ISO 9001 and ISO 13485.



- Equipment and Patient Drapes
- Surgical Products and Accessories
- Anaesthesia and Airway Management
- Veterinary Products

As part of our continuous improvement we also work with Companies to develop new, innovative products, as well as supplying OEM products. Our Contract Manufacturing services include cleanroom moulding, extrusion and product assembly.

CONTENTS

CAMERA DRAPES	04
MICROSCOPE DRAPES	08
ULTRASOUND PROBE DRAPES	12
C-ARM DRAPES	14
OTHER EQUIPMENT DRAPES	17
PATIENT DRAPES	19
SURGICAL ACCESSORIES	23
PROCEDURE PACKS	26

P3 EQUIPMENT DRAPES

All our equipment drapes can also be supplied in bulk, non-sterile.

Latest Catalogues and further product information, including product videos and data sheets can be found on our Distributor Portal; https://p3medical.box.com/v/TrainingVideos

SJM 10/17_1

CAMERA DRAPES



HIGH-DEFINITION CAMERA DRAPE

The next level of camera drapes.

- Suitable for high-definition camera systems without compromising the HD image quality.
- Closed-end offers camera full protection from fluids and acts as a reliable sterile barrier.
- Allows quick and easy scope changes during procedures without compromising sterility.
- Broad-edged ring ensures ease of use and prevents contamination of the sterile field.
- Sterile, single-use.

Size 150mm x 2.5m,	Product code	Unit of sale	
ECN150-250CAP Ring aperture 85mm 50 unit	ECN150-250CAP	50 units	



01 Attach the endoscope Simply push the endoscope into the drape holder.



2 Attach the camera

Connect the camera onto the camera drape holder and extend the drape over the cable by pulling back the ring. Any excess drape is retained inside the ring.

QUICK CHANGE CAMERA DRAPE • Simple and secure 2-step fixing method. TRAINING • Allows quick and easy scope changes during procedures VIDEOS without compromising sterility. AVAILABLE • Closed-end offers cameras full protection from fluids and acts as a reliable sterile barrier. • Broad-edged ring reduces risk of contaminating the sterile field. • Sterile, single-use. **Dimensions Product code** Unit of sale Size 150mm x 2.5m, ECE150-250C 50 units Ring aperture 85mm Attach the camera

Simple push-fit. The drape can then be extended over camera and cable by simply pulling back the ring. Any excess drape is retained inside the ring.



CONCENTRIC RING STYLE

- Ring style drapes ensure easy and efficient application without compromising sterility.
- Drape is fully contained within the ring and can be simply pulled back over camera and cable.
- Any excess material remains securely inside the ring.
- Available with 3 different closures.
- Sterile, single-use.

Product code	Dimensions	Unit of sale
ECN150-250	with two adhesive tape closures Size 150mm x 2.5m, Ring aperture 85mm	50 units
ECN150-250S	elastic collar with two twist-lock spring clips and one adhesive tape closure Size 150mm x 2.5m, Ring aperture 85mm	50 units
ECN150-250TK	elastomeric end section with two adhesive tape closures Size 150mm x 2.5m, Ring aperture 85mm	50 units





TELESCOPIC STYLE

- Folded cuff to facilitate application without compromising sterility.
- Sterile, single-use.

Product code	Dimensions	Unit of sale
E0D130-250T	with perforated tip and two adhesive tape closures Size 130mm x 2.5m	25 units
E0D165-250TK	with elastomeric end section and one adhesive tape closure Size 165mm x 2.5m, Aperture 10mm dia	25 units
E0D130-250TB	with cardboard sleeve Size 130mm x 2.5m	25 units



P³

CONCERTINA STYLE

- Locator insert assists hand access through the drape.
- Folded cuff to facilitate application without compromising sterility.
- Sterile, single-use.

Product code	Dimensions	Unit of sale
E0D150-300	supplied with two adhesive tape closures Size 150mm x 3m	25 units
EOD165-300	supplied with two adhesive tape closures Size 165mm x 3m	25 units

FLAT FOLDED STYLE

- Folded cuff to facilitate application without compromising sterility.
- Sterile, single-use.

Product code	Dimensions	Unit of sale
E0D130-250	supplied with one adhesive tap closure Size 130mm x 2.5m	e 25 units



DRILL DRAPE – TELESCOPIC STYLE

- Folded cuff to facilitate application without compromising sterility.
- Sterile, single-use.

Product code	Dimensions	Unit of sale
E0D130-100T	supplied with two adhesive tape closures Size 130mm x 1m	20 units

LAPAROSCOPIC INSTRUMENT POUCH


- Specially designed sterile pouch with five pockets to provide temporary storage for instruments during surgical procedures.
- Sterile, single-use.

Product code	Dimensions	Unit of sale
LIP500-048	with adhesive area to attach pouch to drapes/surfaces Size 500 x 480mm	25 units

TRAINING VIDEOS AVAILABLE

MICROSCOPE DRAPES





STANDARD

MICROSCOPE DRAPES

- To fit specific objective lenses.
- With Velcro securing tapes at each eyepiece and two/three gathering straps for a neat and secure fit.
- With removable objective lens cover.
- Sterile, single-use.

To fit ZEISS OPMI microscopes (objective lens 48mm o/dia)

Product code	Dimensions	Unit of sale
M03110-160	Size 1.35 x 1.6m (standard)	10 units
M03110-200	Size 1.35 x 2m (long)	10 units
M03110-260	Size 1.35 x 2.6m (extra-long)	10 units


To fit ZEISS MD microscopes (objective lens 65mm o/dia)

Product code	Dimensions	Unit of sale
MD3110-160	Size 1.35 x 1.6m (standard)	10 units
MD3110-200	Size 1.35 x 2m (long)	10 units
MD3110-260	Size 1.35 x 2.6m (extra-long)	10 units


To fit LEICA/WILD microscopes (objective lens 70mm o/dia)

Product code	Dimensions	Unit of sale
MW3110-160	Size 1.35 x 1.6m (standard)	10 units
MW3110-200	Size 1.35 x 2m (long)	10 units
MW3135-240	Size 1.35 x 2.6m (extra-long)	10 units


TRAINING VIDEOS AVAILABLE



LONG



EXTRA LONG



MICROSCOPE DRAPES, UNIVERSAL FITTING

- With elastomeric panel (140mm o/dia, 40mm i/dia) to fit more than one objective lens size.
- Velcro securing tapes at each eyepiece and two/three gathering straps for a neat and secure fit.
- Sterile, single-use.

Product code	Dimensions	Unit of sale
MC3106-160	Size 1.35 x 1.6m (standard)	10 units
MC3106-260	Size 1.35 x 2.6m (extra-long)	10 units



EXTRA LONG

10

OBJECTIVE LENS COVER

- With removable lens protector.
- For use in conjunction with the P3 universal fitting microscope drape.
- Sterile, single-use.

Product code	Dimensions	Unit of sale
MLC 048	48mm dia (for Zeiss OPMI)	10 units
MLC 065	65mm dia (for Zeiss MD)	10 units
MLC 070	70mm dia (for LEICA /WILD)	10 units




P³

MICROSCOPE HANDLE DRAPES

- Clear plastic with adhesive fixing tape.
- Packed in pairs.
- Sterile, single-use.

Product code	Dimensions	Unit of sale
MHD180-330	Size 180 x 330mm	50 pairs





STERILE Aquasonic

RANSANSSION GEL WATER SOLUBLE

PARKER LABORATORIES INC

NTEED STERILE UNLESS

CE 0086





ULT127-150T



ULTRASOUND PROBE DRAPES

- Reduces the risk of patient-to-patient transmission of infectious agents.
- Sterile, single-use.

Telescopic style

Product code	Dimensions	Unit of sale
ULT127-150T	with two adhesive tape closures Size 127mm x 1.5m	25 units
ULT127-150T-G	supplied with sterile gel sachets Size 127mm x 1.5m	25 units


TRAINING VIDEOS AVAILABLE

ULT127-150T-G

Flat folded style

UGD130-082 with two adhesive ta Size 130 x 820mm	pe closures 25 units





C-ARM FLUOROSCOPE DRAPE


- Elasticated collar (contains latex) for easy application and increased security.
- Large size to accommodate most C-Arms.
- Sterile, single-use.


Product code	Dimensions	Unit of sale
CFD760-090	Size 760 x 900mm	25 units

C-ARM DRAPE, EXTRA LARGE

- Extra-large size to accommodate mobile C-Arms.
- Supplied with six adhesive tape closures for a secure fit.
- Sterile, single-use.

Product code	Dimensions	Unit of sale
CAD120-150	Size 1.2 x 1.5m	10 units





C-ARM DRAPE, SIDE OPENING


- Side opening to facilitate easy draping.
- Folded cuff to facilitate application without compromising sterility.
- Full-width adhesive tape closure for a secure fit.
- Sterile, single-use.

Product code	Dimensions	Unit of sale
CAD600-100	Size 600mm x 1m	10 units

C-ARM DRAPE, FULL SIZE

- Telescopic style and easily applied by one person.
- Supplied with ten adhesive tape fasteners and four gathering strips for a secure fit.
- Sterile, single-use.

Product code	Dimensions	Unit of sale
CAD265-152	Size 2.65 x 1.52m	10 units



TRAINING VIDEOS AVAILABLE




BONNET/BANDED BAG

- Elasticated collar for easy application and increased security.
- With many uses in operating theatres and X-ray.
- Sterile, single-use.

Product code	Dimensions	Unit of sale
CAB070	700mm dia	25 units
CAB080	800mm dia	25 units
CAB090	900mm dia	25 units
CAB100	1.0m dia	25 units







X-RAY SCREEN COVER


- Transparent film to allow clear view of screen during sterile procedures.
- Supplied with four adhesive tapes (two full-width) for a secure fit.
- Sterile, single-use.


Product code	Dimensions	Unit of sale
SCC090-120	Size 900 x 1200mm	25 units

MEMBRANE KEYPAD DRAPE

- Simple clear cover for enhanced security when using equipment keypad during sterile procedures.
- Supplied with four adhesive fixing tapes.
- Sterile, single-use.

Product code	Dimensions	Unit of sale
KPD200-210	Size 200 x 210mm	25 units







COLPOSCOPE DRAPE, UNIVERSAL FITTING

- Clear plastic with elastomeric panel (140mm o/dia, 40mm i/dia) to fit more than one objective lens size.
- Single-use, packed under cleanroom conditions

Product code	Dimensions	Unit of sale
ID0001NS	Size 840 x 760mm	100 units



•З



UROLOGY DRAPE


- Blue plastic film for excellent shape conformity.
- 50mm wide adhesive strip with finger lift for easy application and positioning.
- Latex finger cot with additional adhesive ring around the base.
- Aperture 60mm dia.
- Sterile, single-use.


Product code	Dimensions	Unit of sale
UR0600-090	Size 600 x 900mm	10 units
UR0600-126	Size 600mm x 1.26m	10 units

U-DRAPE

- High-quality blue plastic film.
- Adhesive tape with finger lift for easy application.
- Sterile, single-use.

Product code	Dimensions	Unit of sale
UAP130-120	Size 1.2 x 1.3m	25 units





LOWER LIMB EXCLUSION DRAPE


- High-quality clear plastic film.
- Folded cuff has full width 50mm adhesive strips with finger lift for easy closure.
- Sterile, single-use.


Product code	Dimensions	Unit of sale
LLD400-700	Size 400 x 700mm	10 units

MAXILLO-FACIAL DRAPE

- High-quality blue plastic film for excellent shape conformity.
- With latex finger cot.
- Full width 50mm adhesive tapes with finger lift and adhesive ring surrounding each aperture.
- Sterile, single-use.

Product code	Dimensions	Unit of sale
MXF600-950	Size 600 x 950mm	10 units





SMALL APERTURE DRAPE

- High-quality blue plastic film for excellent shape conformity.
- 100mm square adhesive tape around aperture.
- With finger lift for easy application.
- Sterile, single-use.


21

Product code	Dimensions	Unit of sale
SMA400-040	Size 400 x 400mm, Aperture 60mm dia	25 units

LARGE DRAPE WITH APERTURE AND POUCH

- High-quality blue plastic film with fluid collection pouch.
- 156mm square adhesive area around aperture.
- Sterile, single-use.


Product code	Dimensions	Unit of sale
MAS100-220	Size 1.0 x 2.2m, Aperture 60mm dia	10 units
MAS100-220	,	10 units



TOURNIQUET DRAPE

- High-quality blue plastic film.
- Suitable for both upper and lower limb procedures.
- With two full width adhesive strips (25mm and 30mm) for maximum effectiveness.
- Designed specifically for tourniquet cuff based procedures to keep the padding and tourniquet cuff clean and dry during pre-op skin preparation.
- Single-use, packed under cleanroom conditions.

Product code	Dimensions	Unit of sale
TRN090-300NS	Size 300 x 860mm	100 units





ELASTOMERIC TOURNIQUET DRAPES

- High-quality clear plastic film with elastomeric panel for quick and efficient application.
- Suitable for both upper and lower limb procedures.
- Single-use, packed under cleanroom conditions.

Product code	Dimensions	Unit of sale
TRA060NS	Size 550 x 650mm, Aperture 60mm dia	50 units
TRA075NS	Size 550 x 650mm, Aperture 75mm dia	50 units
TRA100NS	Size 550 x 650mm, Aperture 100mm dia	50 units

ABSORBENT PROCEDURE DRAPE

- High-quality blue plastic film with white absorbent pad to absorb up to 1 litre of fluid.
- Drier floors lead to a safer environment and faster patient turnaround.
- Sterile, single-use.

Product code	Dimensions	Unit of sale
APD090-125	Size 900mm x 1.45m, Aperture 120 x 60mm (150mm using perforation) Absorbent pad size 560 x 820mm	25 units






DISPOSAFE SHARPS MAT

- Safe and cost-effective method for sharps disposal.
- Size (before folding) 145 x 120mm.
- Sterile, single-use, latex free.

Product code	Unit of sale
DS-599CTN	50 units (one dispenser carton)
DS-599	300 units (one case of 6 dispenser cartons)



DISPOSAFE PLUS

- Sharps mat with enhanced safety features.
- Tough plastic backing gives increased resistance to penetration.
- Extra-large adhesive area gives enhanced visibility and holds more sharps, with separate adhesive strip for blades.
- Size (before folding) 145 x 120mm.
- Sterile, single-use, latex free.

Product code	Unit of sale
QF-1001CTN	50 units (one dispenser carton)
QF-1001	300 units (one case of 6 dispenser cartons)



CAUTERY TIP CLEANING PAD

- Tip cleaning pad with X-Ray detectable strip.
- Full width adhesive tape on reverse, with finger-lift for secure positioning.
- Size 50 x 50mm.
- Sterile, single-use, latex free.

Product code	Unit of sale
CTC050	50 units






SURGICAL SKIN MARKER

- Non-toxic, prep resistant gentian violet ink can be removed with isopropyl alcohol.
- Versatile tip for thin or broad line marking.
- Each pen supplied with flexible ruler (centimetres).
- Sterile, single-use, latex free.

Product code	Unit of sale
CSMK-CTN	Box of 10 units
CSMK	300 units (30 boxes of 10 units)





LAPAROSCOPIC PROCEDURE PACK

- Designed for laparoscopic procedures.
- Reduces preparation time.
- Saves storage space.
- Contents can be customised on request.
- Supplied sterile, single-use.

Set can include insufflation tubing, choice of P3 camera drapes, sharps pad or box, swabs, blade, etc.

Details of options available on request.

UROLOGICAL PROCEDURE PACK

Designed for use in the following procedures:

- Flexible cystoscopy
- Bladder pressure studies
- Supra pubic catheter insertion
- Urethral dilatation examination
- Prostate microwaving

Set includes large fenestrated drape with 1 litre capacity absorbent pad, plus customised selection of other equipment to end-user specification.

Items can include gauze swabs, gallipots, kidney dish, paper drapes, forceps, irrigation tubing, etc.

Supplied sterile, single-use. Details of further options available on request.





HSG PROCEDURE PACK

Designed for use in Hysterosalpingogram procedures.

- Designed for laparoscopic procedures.
- Reduces preparation time.
- Saves storage space.
- Contents can be customised on request.
- Supplied sterile, single-use.

Standard set includes patient drape with fenestration and absorbent pad, green crepe outer, 250ml gallipot, 2 Rampley forceps and 5 gauze swabs.

Vaginal speculum and kidney dish can be included as optional extras.

Details of further options available on request.



P3 Medical Limited 1 Newbridge Close Bristol BS4 4AX United Kingdom

T: +44 (0) 117 972 8888 F: +44 (0) 117 440 9493 E: export@p3medical.com